


PyOptomip User's Guide

12/12/2022



Purpose

PyOptomip is a python application which aims to facilitate testing of silicon electronic photonic integrated circuits by controlling testing equipment through a user interface.

Table of Contents

1. Connecting Instruments

Getting Started

IDA Stage Instrument Connections

2. Manual Measurements

Instrument Control

Home Tab

Electrical Tab

Optical Tab

3. Automated Measurements

Labelling Devices

Generating Coordinate Files

Creating Test Routines

Importing Testing Parameters

Filtering Devices

Creating Transform Matrices

Setting Wedge Probe Limit

Starting Measurements


4. Advanced

Requirements

PyOptomip Github Repo

Getting Started

To begin using PyOptomip to perform measurements you'll first need to connect any instruments you will be using.



1. Select the device you wish to connect to from the drop-down menu on the connect instrument page.
2. Once the desired device is selected, enter its GPIB address in the designated space.
3. After all necessary instruments are connected, press done.

IDA Stage Instrument Connections

The GPIB addresses for the instruments at the IDA stage in KAISER 4060 are as follows:

Stage: Corvus Eco - ASRL7::INSTR


Stage: Electrical Probing - ASRL10::INSTR

Laser: hp816x with N77 Detector -

GPIB0::20::INSTR

USB0::0x0957::0x3718::MY48102149::INSTR

Source Meter Unit - GPIB0::26::INSTR



Instrument Control

PyOptomip has five tabs which are used to control instruments and automate measurements.

Home tab:

Stage control, detector measurements and settings

Electrical tab:

SMU control, IV sweeps

Optical tab:


Laser control, wavelength sweeps

Automated Measurements:

Automated measurement settings


Testing Parameters:

Testing routine creation



Home Tab

The home tab can be used to move the stages on which the chip and wedge probe rest, respectively. It can also be used to gauge alignment with a device using the detector and SMU as well as the camera.






Camera Settings

The Camera Settings can be used to adjust the exposure and saturation of the camera and to switch between views.

Begin by pressing **Open** and then use the **Switch Views** button to toggle between a top and side view of the chip.

Make sure to close the camera before exiting PyOptomip.

If you wish to perform a video of your measurements simply choose a folder in which to save your results and then use the **Start** and **Stop** buttons.



Corvus Eco

The three axes of the stage on which the chip rests can be controlled by entering a number of micrometres and pressing either the + or - button for that axis.

Select **Auto measurement** for the detector channels whose readings you want to see. You can use these readings to determine when the laser is more or less aligned.

At this point, **Fine Align** can be used to optimize the light coupled to the detector.

The screenshot displays the Corvus Eco control interface. At the top, under the 'Corvus' header, there are three 'Motor Control' rows. Each row has a minus button, a text input field containing the value '1', '2', and '3' respectively, followed by 'um' and a plus button. Below this is a 'Fine align' section with two tabs: 'Fine Align' (selected) and 'Settings'. Under the 'Fine Align' tab, there is a 'Detector Settings' section. It includes two input fields: 'Initial range (dBm)' with the value '-20' and 'Range dec. (dBm)' with the value '20'. Below these are four rows for detector settings. Each row has a checkbox for 'Auto measurement' (all are unchecked), a checkbox for 'Include in sweep' (all are unchecked), and a 'Power (dBm):' field with the value '-100'. The rows are labeled 'Slot 1 Det 1', 'Slot 2 Det 1', 'Slot 3 Det 1', and 'Slot 4 Det 1'.

Corvus		
Motor Control	-	1 um +
Motor Control	-	2 um +
Motor Control	-	3 um +
Fine align		
Fine Align Settings		
Detector Settings		
Initial range (dBm)	-20	Range dec. (dBm) 20
Slot 1 Det 1	<input type="checkbox"/> Auto measurement	<input type="checkbox"/> Include in sweep Power (dBm):-100
Slot 2 Det 1	<input type="checkbox"/> Auto measurement	<input type="checkbox"/> Include in sweep Power (dBm):-100
Slot 3 Det 1	<input type="checkbox"/> Auto measurement	<input type="checkbox"/> Include in sweep Power (dBm):-100
Slot 4 Det 1	<input type="checkbox"/> Auto measurement	<input type="checkbox"/> Include in sweep Power (dBm):-100

Electrical Stage


The stage on which the wedge probe sits can be controlled by entering a number of micrometres and pressing either the + or - button for that axis. The resistance measured by the SMU is displayed to aid in aligning the probe.

Wedge Probe Stage				
X	<input type="button" value="-"/>	<input type="text" value="500"/>	um	<input type="button" value="+"/>
Y	<input type="button" value="-"/>	<input type="text" value="10"/>	um	<input type="button" value="+"/>
Z	<input type="button" value="-"/>	<input type="text" value="100"/>	um	<input type="button" value="+"/>

Electrical Connection Detector			
Resistance (Ω):			
Channel A:	70250000.0	Channel B:	18.5032

Electrical Tab

The electrical tab can be used to control an SMU to set the voltage and current on its channels. It can also be used to perform IV sweeps and plot the results.



SMU Control

Channel A, or B, or both, of the SMU can be set by checking off the desired channel, inputting parameters and pressing **Set**.

Once parameters are set, the selected channels can be turned on.

Readings from each channel are displayed at the bottom of the panel.

SMU Control			
Select SMU Output	<input checked="" type="checkbox"/> A	<input checked="" type="checkbox"/> B	ON/OFF
Set Voltage (V)	0		Set
Set Current (mA)	0		Set
Set Voltage limit (V)	0		Set
Set Current limit (mA)	0		Set
Set Power limit (mW)	0		Set
	Channel A	Channel B	
Voltage Reading (V):	-5.8651e-06	0.000953674	
Current Reading (mA):	-2.2068e-10	-1.81198e-12	
Resistance (Ω):	26577.397	-526315963.752	


Sweep Settings

IV Sweeps can be performed with channel A or B by checking their respective boxes and clicking **Set**.

Choose the minimum and maximum Voltage in *Volts* as well as the resolution in *Volts*.

If you wish to save the sweep results You can choose a folder to do so.

You can also choose which channel to plot as well as the type of plot. When ready, press the **IV Sweep** button.




The image shows a software window titled "Sweep Settings". It contains several input fields and checkboxes. At the top, "Select SMU Output" has two checkboxes, "A" and "B", both of which are checked. To the right of these is a "Set" button. Below this are three input fields: "Set Max. Voltage (V)" with the value "7", "Set Min. Voltage (V)" with the value "1", and "Set Resolution" with the value "12". Underneath these is a "Save folder:" label followed by a text box containing "D:\\Work\\PyOptomip_testdata" and an "Open..." button. The "Plot Selection:" section has checkboxes for "A" (checked), "B" (unchecked), "IV" (unchecked), "RV" (unchecked), and "PV" (checked). At the bottom of the window is a large button labeled "IV Sweep".

Parameter	Value / Selection
Select SMU Output	<input checked="" type="checkbox"/> A, <input checked="" type="checkbox"/> B
Set Max. Voltage (V)	7
Set Min. Voltage (V)	1
Set Resolution	12
Save folder:	D:\\Work\\PyOptomip_testdata
Plot Selection	<input checked="" type="checkbox"/> A, <input type="checkbox"/> B, <input type="checkbox"/> IV, <input type="checkbox"/> RV, <input checked="" type="checkbox"/> PV
IV Sweep	[Button]

Optical Tab

The optical tab can be used to control the laser as well as perform and plot optical spectrum sweeps.



CW Settings

The laser output parameters can be set using the CW settings

The laser slot should be set to auto by default and should in most cases stay this way.

Use the laser status **ON OFF** buttons to toggle the laser On and Off

Select the desired laser output and that click the **Set** button to push output to laser.

The power and wavelength of the laser can be set in dBm and nm respectively.


CW Settings		
Laser slot:	auto ▾	
Laser status:	OFF	ON
Laser output	High power ▾	Set
Power (dBm)	0	Set
Wavelength (nm)	0	Set

Sweep Settings

In the sweep settings frame all parameters required for an optical sweep can be set.

Starting wavelength, stop wavelength and step size between wavelengths can all be set here in nm.

Choose which detectors you wish to include in the sweep in the Home tab.



The image shows a software window titled "Sweep Settings". It contains several input fields and dropdown menus for configuring an optical sweep. The parameters are: Start (nm) set to 0, Stop (nm) set to 0, Step (nm) set to 0, Sweep power (dBm) set to 0, Sweep speed set to "auto", Laser output set to "High power", and Number of scans set to 1. Each of these fields has a small downward arrow icon on its right side, indicating they are dropdown menus. At the bottom right of the window is a button labeled "Sweep".

Parameter	Value
Start (nm)	0
Stop (nm)	0
Step (nm)	0
Sweep power (dBm)	0
Sweep speed	auto
Laser output	High power
Number of scans	1

Sweep


Labelling Devices

In order to perform automated measurements, you'll need to generate a coordinate text file from the GDS of the chip you wish to measure. You'll need to install the [SiEPIC-Tools package](#) for KLayout and the [SiEPIC-EBeam-PDK](#) to do this.

For every device in your design ensure you place labels for the optical and electrical connections as necessary. Make sure to follow the labelling standards and design standards outlined in the next pages.

Optical Labels


Optical labels should be placed at the tip of the grating coupler to be connected to the laser. Labels should follow the format ***opt_<polarization>_<wavelength>_<device>_<unique-device-id>***. Each electro-optic device should have an id which is unique from any other on the same chip.



opt_@opt_TE_1550_device_<unique_device_id>


Electrical Labels

Electrical labels should be placed at the center of each bond pad and should follow the format ***elec_<unique-device-id>_<pad-name>***.




Electrical Pad Arrangement

Bond pads should be $75\mu\text{m}$ by $75\mu\text{m}$ and be spaced $125\mu\text{m}$ apart. There should be at least two bond pads and the leftmost pad should be signal ground. Ensure that all bond pads for a device are above the associated grating couplers in the orientation shown here and that there is adequate space between grating couplers and bond pads.




Generating Coordinate Files

Once you've finished your design and labelled all devices you can export a text file of device coordinates by clicking on the **SiEPIC** toolbar and then selecting **Layout** followed by **Automated Coordinate Extraction**. A pop up window will open which will allow you to review and then save your text file.



Creating Test Routines

Before the user can begin running automated tests they must create a set of routines for each device that they want to test. This is done using the testing parameters tab in pyOptomip or by running the standalone testing parameters python script. The left half of the tab can be used to create routines while the right half can be used to associate routines with devices.



Uploading Device Coordinate Files

Instrument Control

Home Electrical Optical Automated Measurements Testing Parameters

Test Parameter Creation

No file selected Open...

Routine Creation

Routine Subroutine

Routine Name:

Set Max (V/mA):

Set Min (V/mA):

Set Resolution (mV/mA):

Plot Type: ☐ IV/VI ☐ RV/RI ☐ PV/PI

Select SMU Channel: ☐ A ☐ B

Start (nm)

Stop (nm)

Step (nm)

Sweep power (dBm)

Initial Range (dBm)

Range Decrement (dBm)

Sweep speed

Laser output

Number of scans

Wavelengths (nm)

Voltages (V)

Save as Subroutine

Device Linking

Select All Unselect All Select keyword Unselect keyword

Select Wavelength: ☒ 1310 ☒ 1550

Select Polarization: ☒ TE ☒ TM

Group Devices Device Data

Associated Routines

Remove Routine

Set

Save folder: Open...

Send to Automeasure Import Export

Use this to open folder and choose chip coordinate file

Creating Routines

Begin by selecting the type of routine you wish to perform. Next, give a name to your routine and populate the parameters of the routine in the given text boxes. Finally, press **Save as Subroutine** and your routine will appear in the Subroutine column of the corresponding routine type.

Routine Creation

Routine

☒ Wavelength Sweep
☐ Voltage Sweep
☐ Current Sweep
☐ Set Wavelength Voltage Sweep
☐ Set Wavelength Current Sweep
☐ Set Voltage Wavelength Sweep

Subroutine

☐ Default

Routine Name:

Start (nm)

Stop (nm)

Step (nm)

Sweep power (dBm)

Initial Range (dBm)

Range Decrement (dBm)

Sweep speed

Laser output

Number of scans


auto

High power

1

Save as Subroutine

Setting Routines



The screenshot shows a 'Device Linking' window with the following components:

- Buttons:** 'Select All', 'Unselect All', 'Select keyword', and 'Unselect keyword'.
- Filters:**
 - 'Select Wavelength:' with checkboxes for 1310 and 1550 (both checked).
 - 'Select Polarization:' with checkboxes for TE and TM (both checked).
- Table:** A table with four columns: 'Group', 'Devices', 'Device Data', and 'Associated Routines'.


Group	Devices	Device Data	Associated Routines
<input type="checkbox"/> device	<input type="checkbox"/> device4 <input type="checkbox"/> device3 <input type="checkbox"/> device2 <input type="checkbox"/> device1		
- Footer:** A 'Set' button on the left and a 'Remove Routine' button on the right.

Next, select routines by checking them off and select the devices you'd like to associate with the routines. You can filter the devices by wavelength, polarization or by selecting the device type in the Group column. Click **Set** to link routines to devices. Once a device has associated routines you can view them by selecting the individual device. If you wish to remove a routine from a device you can use the **Remove Routine** button.

Exporting Routines

If you wish to save your devices and associated routines to perform measurements later, or just to have a back-up file, you can choose the folder in which to save your file and then press **Export**.

If you wish to perform measurements right away, you can send the routines and devices straight to the Automated Measurements tab with **Send to AutoMeasure**.




Save folder: Open...

Send to Automeasure Import Export


Performing Automated Measurements

To automate measurements of devices with the routines you've set in the Testing Parameters tab you'll need to go to the Automeasure tab. Here you can choose which devices to test and begin your automated measurements.



Importing Testing Parameters

If you created your testing parameters file in a separate instance of PyOptomip than the one currently running, you can import your file using the **Import Testing Parameters** button. All the devices should then populate in the device list at the top of the panel. Devices with associated routines will appear in black, while those with no routines will appear in grey. You can then select the devices you wish to test using the **Select All**, **Unselect All**, or **Filter** buttons.




Filtering Devices

Within the filter frame you can select which modes and wavelengths you wish to measure. For example, selecting **TE** will select all TE devices.

You can also select the device types that you wish to include by choosing them from the drop down menu then pressing **Select**.

Keywords can be used to select or unselect any device which contains the desired keyword in its ID.



The screenshot shows a 'Filter' dialog box with the following sections:

- Mode:** Contains two checkboxes, ☐ TE and ☐ TM.
- Wavelength:** Contains two checkboxes, ☐ 1310nm and ☐ 1550nm.
- Device Type:** A dropdown menu with a downward arrow.
- Keywords:** A text input field.
- Buttons:** 'Select' and 'Unselect' buttons are located to the right of the 'Device Type' and 'Keywords' fields respectively. A 'Done' button is at the bottom center.

Creating Transform Matrices

Before automated measurements can commence the transform matrices must be created. Within the automated measurements tab are the electrical and optical alignment boxes. Align the fibre array manually (through the home tab) onto a device. Once the array is aligned, choose the device you have aligned to from the drop down menu and hit the **Get Pos.** button to save the current coordinates. Repeat this until all boxes have been filled. For the optical matrix there is the option to save the alignment and import it later.

Choose alignment device
from dropdown

		X	Y	Z	
Device 1	device1	4537.55188	4480.51453	-744.62457	Get Pos.
Device 2	device2	4592.55981	4545.53223	-744.62457	Get Pos.
Device 3	device3	4902.55737	3980.52979	-744.62457	Get Pos.
Save Optical Alignment		Import Optical Alignment			

Once the array is aligned click the
Get Pos. Button

Setting Wedge Probe Limit and Scaling

Move the wedge probe as close as you're comfortable getting it to the fibre array and then press the **Set Position** button to ensure the wedge probe never moves past this position.

If the axes of the wedge probe and stage are slightly misaligned you may need to adjust the x and y scales.


Minimum Wedge Probe Position in X:			Scale Adjust		
<input type="text" value="-250"/>	<input type="button" value="Set Position"/>	<input type="button" value="Undo"/>	X: <input type="text"/>	Y: <input type="text"/>	<input type="button" value="Set Scale Adjustment"/> <input style="border: 2px solid blue;" type="button" value="?"/>

Starting Measurements

After uploading the coordinate file, creating the transform matrices and limiting the wedge probe movement you are ready to begin taking automated measurements.

Select which detectors you wish to use for measurements and then set the folder in which you wish to save your measurement data by clicking **Open...** beside the designated text box.

Finally, select **Start Measurements**.



The screenshot shows a software dialog box titled "Choose Detectors for Automated Measurements". It contains several controls for configuring automated measurements:

- Detectors:** Four checkboxes are present: "Slot 1 Det 1" (checked), "Slot 2 Det 1" (unchecked), "Slot 3 Det 1" (unchecked), and "Slot 4 Det 1" (unchecked).
- Save Results:** A section with the label "Save Results" and "Save Folder:".
- Save Folder:** A text input field containing the path "C:\Users\SiEPIC_Kaiser\Nextcloud\Lab data", followed by an "Open..." button.
- Electro-Optic Measurements:** A section with the label "Electro-Optic Measurements" and "Align to Device".
- Align to Device:** A dropdown menu (currently showing a downward arrow) and a "Go" button.
- Start Measurements:** A large, prominent button at the bottom of the dialog.

PyOptomip Github Repo

To make changes to the PyOptomip code or to view how the software works, all related PyOptomip content can be found in the following github repo:

<https://github.com/SiEPIC/pyOptomip>

The ida-main branch is the version currently running on the IDA stage in KAISER 4060.

Requirements

Before one can begin using PyOptomip on their PC they must ensure their computer has the required modules and API's. To find all required python add-ons and modules use the **requirements.txt** file located in the PyOptomip folder in the Github repository.